



TECHNICAL DATA

Rev. 03.04.2020

CODE	DESCRIPTION	UNITS (Metric)	DSU MOBILE
1.0	COLD WATER NETWORK		
	Built-in and easily accessible shut-off valve	---	•
	Connection	"G	1
	Pressure (min. - max.)	bar	4 - 8
5.0	GULLY DRAIN/FLOOR SYPHONIC DRAINAGE		
	Connection (internal diameter)	mm	80
	Maximum Temperature	°C	100
	Temperature resistant (copper, cast iron or stainless steel)	---	•
	Built with anti-corrosive materials	---	•
6.0	MAIN ELECTRICITY		
	Protected electrical connection according to current regulations.	---	•
	Type	---	3~+ PE
	Voltage	V	400
	Power frequency	Hz	50
	Power	kW	23
WEIGHTS			
	Weight	Kg	3600
NOTES			
	Connections must be installed in accordance with the regulations in force.		•
	IMPORTANT! A perfectly level support surface is required at the destination. A 200 mm thick concrete floor is recommended.		•

Format A3 Scale 1:25 Date 03.04.2020

Drawn A. MERA

Designed A. MERA

Authorized

Customer

Project

DSU MOBILE

Layout

DISTRIBUCION

Project num.

M-0003-20

Order num.



ANTONIO MATACHANA, S.A. claims intellectual property of this plan/project and reserves the exclusive right to use it. No part of it may be used or copied without the prior permission of the legal representative of Antonio Matachana